

EzQuant ssDNA Quantification Kit

Product Description

EzQuant ssDNA Quantitative Kit is for ssDNA or oligonucleotides. It enables quick and accurate quantification of low-concentration samples from 50 pg/ μ l to 200 ng/ μ l when using with EzCube fluorometer.

It does not specifically recognize ssDNA but can also bind to dsDNA and RNA. However, assay results will not be affected by common contaminants; such as proteins, salts, and detergents.

Components

Material	Concentration	Amount T00-FPBR01-00 (100 Assays)	
EzQuant ssDNA Buffer	Not applicable	50 ml *1	
EzQuant ssDNA Reagent	200X in DMSO	250 µl *1	
EzQuant ssDNA Standard #1	0 ng/µl in TE Buffer	1 ml *1	
EzQuant ssDNA	20 ng/µl	1 ml*1	
Standard #2	in TE Buffer		

Note: For Research Use Only. Not for use in diagnostic procedures.

Shipping and Storage

All the components are shipped with ice pack and can be stored at 2 - 8 °C away from light for one year.

Required Material (not provided)

- 1. P2/P20/P200 Micropipette
- 2. Tips
- 3. 0.5 ml thin wall polypropylene tubes
- 4. 1.5 2 ml microcentrifuge tubes or 15/50 ml centrifuge tubes

Perform Assay Procedure

- Experimental preparation
- 1. Equilibrate the material to the room temperature before use.
- Prepare enough 0.5 ml thin wall polypropylene tubes and label the tube lids.
 EzQuant ssDNA Assay requires 2 standards.
 Note: Please do not label the side wall of the tube as this could interfere with



the collection of fluorescent signals.

Fresh prepare the 1X Working Solution by diluting the EzQuant ssDNA Reagent
 1:200 in EzQuant ssDNA Buffer.

For example, 2 standards and 5 samples, each tube requires about 200 μ l of working solution, so a total of ~1400 μ l is required. (7 μ l of EzQuant ssDNA Reagent + 1393 μ l EzQuant ssDNA Buffer)

Prepare standards

- 1. Add **190 µl** of **1X Working Solution** to the appropriate tube.
- 2. Add **10 μl** of **EzQuant ssDNA Standard #1** and **EzQuant ssDNA Standard #2** to the appropriate tube.
- 3. Vortex gently for 2-3 seconds and avoid creating air bubbles.

Note: If you do not need to establish a standard curve, you can skip this step.

• Prepare samples

- 1. Add **180 199 μl** of **1X Working Solution** to the appropriate tube.
- 2. Add **1 20 µl** of **samples** to the appropriate tube such that the final volume is 200 µl.
- 3. Vortex gently for 2-3 seconds and avoid creating air bubbles.

• Detection

- 1. Incubate at room temperature for 2 minutes in the dark.
- 2. Follow procedure appropriate for your instrument.




Figure 1. Workflow using EzQuant ssDNA Quantification Kit.



Contaminant	Concentration in 10 µl sample	Detection Results
Ammonium sulfate	20 mM	OK
Sodium chloride	200 mM	OK
Ethanol	10%	OK
Triton X-100	0.1%	OK
Bovine serum albumin	1 mg/ml	OK
dNTPs	2 mM	OK
EDTA	2 mM	OK

Table 1. Effects of contaminants on the results of EzQuant ssDNA Quantification Kit

Ordering Information

Catalog No.	Description
BRFP-0200	EzCube Fluorometer (Blue & Red)
BRFP-0300	EzCube Fluorometer (Blue, Red & Green)
T00-FPBR00-00	EzQuant 1X dsDNA HS Quantification Kit (100 Assays)
T00-FPBR04-00	EzQuant 1X dsDNA HS Quantification Kit (500 Assays)
T00-FPBR01-00	EzQuant ssDNA Quantification Kit (100 Assays)
T00-FPBR02-00	EzQuant RNA HS Quantification Kit (100 Assays)
T00-FPBR03-00	0.5 ml Thin Wall PCR Tube

5F, No. 2, Aly 2, Siwei Ln, Zhongzheng Rd, Xindian Dist, New Taipei City 231, Taiwan (ROC) TEL: +886-2-8219-1206 FAX: +886-2-8219-1205 EMAIL: sales@blue-raybio.com

P/N: 401-FPBR02-10